



Viscolite VL700-S
(Standard Model)

Range: 0 to 10,000 cP
Repeatability: Better than 1%
Temperature: -40 to +150°C

Minimum Sample: 50cc (50ml)
Probe Weight: 700g
Display Weight: 500g
Overall Length: 305mm (std)
Materials: 316ss wetted parts
 Acetal enclosure
Environmental Rating: IP65
Safety Certification: Safe Area
Power Supply: 4 x AA cells

Options: Integral Temperature
 Probe and viscosity correction.
 Highly polished sensing surface.



Viscolite VL700-T15
(Small Sample Model)

Range: 0 to 500, 1000 or 5000 cP
Repeatability: Better than 1%
Temperature: -40 to +150°C
 Integral temperature probe and
 viscosity correction

Minimum Sample: 10cc (10ml)
Probe Weight: 500g
Display Weight: 500g
Overall Length: 215mm (std)
Materials: 316ss wetted parts
 Acetal enclosure
Environmental Rating: IP65
Safety Certification: Safe Area
Power Supply: 4 x AA cells

Options: Highly polished sensing
 surface.